



DRX-4C2



DRX-4C2

Fulcrum is proud to present the new TORBAL, Model DRX-4C2 Class II combination Pharmacy scale and pill counter. The DRX-4C2 is NTEP-certified as a prescription scale; it is equipped with a counting feature which is also NTEP-approved and legal for trade. The scale's counting feature is fast, accurate and easy-to-use.



TORBAL has been the leading supplier of Prescription Balances for over seventy years. This scale combines our long tradition of value and reliability with low cost.



Scale Specifications

Model	DRX-4C
Capacity (Max)	300g
Minimum Load (Min)	0.02g
Reading Unit (d)	0.001g
Verification Unit (e)	0.01g
Tare Range	300g
Accuracy Class	II
Temperature Range	+15C to +30C
Weighing Time	<3s
Pan Dimensions	116mm
Scale Dimensions	235 x 245 x 80mm
Power Supply	Input: 120VDC 60Hz 9.5W Output: 12VDC 500mA
Scale's Net Weight (lb/kg)	8.2 / 3.7
Calibrating Weight	300g



Counting Feature Specifications

MW (Minimum piece Weight)	30 mg (3 e)
MSS (Minimum Sample Size)	10 pieces
Sample Size Selection	10, 30, 60



Standard Features



- **Legal for trade with NTEP-approved counting feature**

- Counting transaction receipt printing
- Compounding and recipe-making mode
- Recipe receipt printing
- Large LCD Display



- **RS-232 Interface Ports**

- Accurate and easy-to-use leveling system
- Used capacity indicator
- Stabilization and Net Weight indicators
- Center of zero indicator
- Automatic zero-setting mechanism (AZSM)
 - Initial zero-setting mechanism (IZSM)
- Sealed front panel
 - Draft ring



Legal for trade with NTEP-approved counting feature



The DRX-4C2 is equipped with an NTEP-approved, legal-for-trade counting feature for prescription filling. The counting feature is easy to use and accurate as it counts up to 90 pills per second. The counting feature is designed in the form of a short wizard which leads the user step-by-step through this simple operation, allowing the user to select a preferred sample size for determining average piece weights and to verify the sample weight. The feature is very user-friendly with the counting process usually taking only a few seconds to complete.

*Note: Handbook 44 Sect 2.20 UR.3.12 Correct Stored Piece Weight states "For prescription scales with a counting feature, the user is responsible for maintaining the correct stored piece weight. This is especially critical when a medicine has been reformulated or comes from different lots." Manufacturers do not supply piece weight information for pills. For these reasons, a factory-stored database is not furnished with the unit, and it is strongly recommend the pharmacist builds and updates the database as new supply bottles are used.



Counting transaction receipt printing

This feature allows the user to print a record of a performed counting transaction. The printed receipt indicates the date and time of the performed transaction as well as the average pill weight and actual count. The printed receipt may be used as a logging tool for the counting feature.

11/23/04 10:49pm

Avg. Pill Wt. 0.58[0] g

Pill Count 68 pc



Compounding and recipe-making mode / Recipe receipt printing

The DRX-4C2 is an excellent tool for compounding and recipe making. The compounding feature will add up the recipe ingredients and print a detailed receipt. While compounding, the scale stores the weights of each ingredient in its memory and allows the user at any time to tally the recipe. A printed receipt indicates the date and time of the performed transaction. It will also list in order the weight of each ingredient and indicate the total of the recipe.

11/23/04 10:50pm

INGREDIENT #1 5.00[0] g net

INGREDIENT #2 10.00[0] g net


INGREDIENT #3 20.00[0] g net

TOTAL: 35.00[0] g net





This precision 300-gram capacity Class II scale, with a resolution of .001g, has a tare capacity equal to its range, and it displays a used capacity indicator whenever weighing. The need for a draft shield is eliminated by the unique design of the pan shield.



The die-cast metal base uses a tripod support system with two adjustable legs, which when combined with a bubble level indicator provides the user with an easy-to-level scale.



Calibration, which in many states must be performed by a State Inspector, is now a simple procedure. The scale leads the operator through the process, whose only requirement is a single 300-gram weight. The scale is also configured with all necessary standard sealing features.



Fulcrum Inc.



Fulcrum Inc.
23 Carol Street
Clifton, NJ, 07014
Tel: 1-866-473-6900
Tel: 973-473-6900
Fax: 973-777-8302
E-mail: Sales@Fulcruminc.net
Web: <http://www.fulcruminc.net>